



**AUTOCUT V (VVG)** Specially designed to use with VTEK - VVTEK Vigeo system - MAGNUM® Bard system



**AUTOCUT P (VVG P)** Specially designed to use with VVTEK Vigeo system - PRO-MAG® ULTRA Argon system



**AUTOCUT B (VVE)** Specially designed to use with PRO-MAG® 2.2 Argon system - BIOPTY® Bard system

- Single patient core biopsy needle designed for use with reusable biopsy system
- Available in several needle gauge sizes and lengths
- Centimeter depth markings assist in needle placement
- Attached removable spacer (VVG / VVGP) maintains the distance between the two needle hubs promoting ease of handling and simplifying targeting
- The coaxial outer needle permits multiple biopsies (US and CT guided)
- Colour coded for quick gauge size identification

AUTOCUT Qty/Box 10pcs			
Add <b>C</b> for Coaxial kit ( VVGC / VVGPC / VVEC )			
Add <b>1</b> for Trocar tip ( VVG1 / VVGP1 )			
CODE	SIZE		
● VVG1210	12Gx10cm		
● VVG1410	VVGP1410	VVE1410	14Gx10cm
● VVG1413	VVGP1413	VVE1413	14Gx13cm
● VVG1415	VVGP1415	VVE1415	14Gx15cm
● VVG1610	VVGP1610	VVE1610	16Gx10cm
● VVG1615	VVGP1615	VVE1615	16Gx15cm
● VVG1620	VVGP1620	VVE1620	16Gx20cm
● VVG1810	VVGP1810	VVE1810	18Gx10cm
● VVG1815	VVGP1815	VVE1815	18Gx15cm
● VVG1820	VVGP1820	VVE1820	18Gx20cm
● VVG1825	VVGP1825	VVE1825	18Gx25cm
● VVG2010	VVGP2010	VVE2010	20Gx10cm
● VVG2015	VVGP2015	VVE2015	20Gx15cm

On request Vigeo is available to provide needles with different diameter and length