



- Lightweight and compact to facilitate one-handed operation during procedures
- Precise, rapid-firing mechanism for an intact core tissue sample
- Specimen notch fully exposed when stylet is advanced for ease of sample removal
- The puncture depth can be adjusted between 15 to 22 mm
- Beveled stylet for easy penetration into the lesion with less trauma to surrounding tissues and organs
- Colour coded for quick gauge size identification
- Echogenic tip for enhanced visualization under ultrasound guidance
- Depth markings etched on cannula for precise positioning
- The Coaxial outer needle permits multiple biopsies (US and CT guided)



One-handed design



Coaxial kit ( VVD0C )

### CT-CORE Qty/Box 10pcs

CODE	SIZE	COAXIAL KIT
● VVD01410	14Gx10cm	VVD0C1410
● VVD01415	14Gx15cm	VVD0C1415
● VVD01610	16Gx10cm	VVD0C1610
● VVD01615	16Gx15cm	VVD0C1615
● VVD01620	16Gx20cm	VVD0C1620
● VVD01810	18Gx10cm	VVD0C1810
● VVD01815	18Gx15cm	VVD0C1815
● VVD01820	18Gx20cm	VVD0C1820
● VVD01825	18Gx25cm	VVD0C1825
● VVD02010	20Gx10cm	VVD0C2010
● VVD02015	20Gx15cm	VVD0C2015

On request Vigeo is available to provide needles with different diameter and length