

RELOCK DUO WOMEN'S HEALTH


Preoperative repositionable breast localization device




In compliance with MDR
(Regulation EU 2017/745)




- Preoperative marking of local lesions
- Retractable hooks with unique flexibility
- Echogenic needle tip for ultrasound use
- Hooks release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Fixation piece for preoperative fixing of the wire
- Twisted wire for enhanced visibility and lower risk of accidental cutting
- Double 5mm loop hook




Half exposure



Full exposure



Wire fixation clamp



Mammography

RELOCK DUO Qty/Box 10pcs

CODE	SIZE
VVNSP1910	19Gx10cm
VVNSP2005	20Gx5cm
VVNSP2007	20Gx7cm
VVNSP2010	20Gx10cm
VVNSP2012	20Gx12cm
VVNSP2015	20Gx15cm

On request Vigeo is available to provide needles with different diameter and length