

RELOCK PREMIUM WOMEN'S HEALTH

Preoperative repositionable breast localization device






In compliance with MDR
(Regulation EU 2017/745)


PATENT PENDING




- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Removable handle for mammography procedure
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- 9mm loop hook



Mammography



Wire fixation clamp
with protection tube



Removable handle

RELOCK PREMIUM Qty/Box 10pcs

| CODE | SIZE |
|------------|----------|
| ● VVRP1905 | 19Gx5cm |
| ● VVRP1907 | 19Gx7cm |
| ● VVRP1909 | 19Gx9cm |
| ● VVRP1912 | 19Gx12cm |

On request Vigeo is available to provide needles
with different diameter and length